

Name: _____ Date: _____ Period: _____ Page:

Hydrometer Lab

Objective:

Density and Buoyancy-

*

*

*

Hydrometer-

*

*

*

Procedure-

1. Fill a large plastic beaker with 400 mL of water.
2. Place the hydrometer in the water. Record the water level on the hydrometer scale.
3. Remove the hydrometer and mix in 10 g of salt into the water. Stir until it is dissolved.
4. Place the hydrometer in the salt water and record the water level on the hydrometer.
5. Repeat steps 3 and 4 four more times, each time adding 10 more grams of salt.

